

MBBS

MORTAR FOR BIO-SOURCED BLOCKS

MORTAR SPECIALLY DESIGNED FOR BIO-SOURCED BLOCK MASONRY (HEMP, FLAX, RAPESEED...)

- > Ready mix Masonry mortar with hydraulic lime. Offering the ease and safety use, the flexibility and durability.



PREPARING THE BEDDING MORTAR

- > Mix one bag of MBBS with 3.5 to 4 litres of water.
- > The mortar can be made either in a concrete mixer or using a whisk mixer (slow speed).
- > Mixing time 3 to 5 minutes.

CONSUMPTION AND PERFORMANCE

	Quantity / m ² of blocks in traditional laying	Quantity / m ² of blocks laid with notched trowel
10 cm block	7 kg / m ²	3,5 kg / m ²
15 cm block	10,5 kg / m ²	5,5 kg / m ²
20 cm block	14 kg / m ²	7 kg / m ²
30 cm block	21 kg / m ²	11 kg / m ²

SHELF LIFE & GUARANTEE

One year from production date if protected in the original packaging and stored in dry conditions. Manufacturer civil Responsibility.

WORKING TEMPERATURE

Use between 5°C and 30°C.

ADDITIONAL INFORMATION

- > Water retention > 90%
- > Setting time : > 3 hours

PACKAGING

- > 25 kg bags. 40 bags per pallet (1T).



by

SAINT-ASTIER
LIME, LIFELONG EXCELLENCE

Chaux de Saint-Astier
28 bis route de Montanceix - La Jarthe - 24110 Saint-Astier
www.saint-astier.co.uk - contact_uk@saint-astier.com
Follow us on @ f in v

